

BRYANRADIOLOGYASSOCIATES

Brazos Valley's leader in advanced diagnostic imaging since 1962
2722 Osler Blvd Bryan, TX 77802; 979-776-8291

Breast & Body MR PROTOCOLS

Breast

Ax Loc
3 Plane Loc
Ax Shim
Map
Ax T1
Ax T2 Fsat
3D pre
3D Dynamic
Sag-High Resolution-Affected breast
Ax T2-High Resolution-Bilateral

Abdomen

3 Plane Loc
Shim
Map
Ax T2 BH
Ax T2 Fs BH
Ax Out phase
Ax In phase
Ax In phase Fsat
Cor SSFP BH
Cor T2 Fase BH (Haste)

Abdomen (Liver)

3 Plane Loc
Map
Shim
Cor T2 FASE
Ax FASE
Ax Out phase Fs
Ax In phase Fs
Cor 3D-Post
Q3D-Ax T1 Fs-Post

Abdomen (Renal Mass)

3 Plane Loc
Shim
Map
Cor T2 FASE BH
Ax T2 FASE
Ax Out phase
Ax In phase
Cor 3D-post
Q3D Ax T1 Fs-post

Adrenals

Loc 3 AX
Shim
Map
Cor T2 FASE Fsath BH
Ax T2 BH
AX T2 Fs BH
Ax In phase BH
Ax Out phase BH
Ax T1 Fs-post
Cor T1 Fs-post
Cor 3D BH post

Uterus

3 Plane Loc
Shim
Map
Sag T2
Ax T2
Cor T2
Ax Out phase
Ax In phase
Ax T1 Fsath-post
Sag T1 post