

BRYAN RADIOLOGY ASSOCIATES

2722 Osler Boulevard Bryan, TX 77802 Phone: (979) 776-8291 Fax: (979) 774-7871

Internal Auditory Canal (IACS)

Indications :

Non Contrast: Cholesteatoma, Hearing Loss, Mastoiditis

| | |
|--------------------------|---|
| Scan Type | Helical |
| Scan Start\ End Location | Scan 1cm below base of skull through mastoid air cells |
| Scan Details | 1mm X 1mm in bone high resolution and soft tissue windows |
| IV Contrast | N/A |
| Scan Delay | N/A |
| Reformat Details | Coronals, magnified views of right and left IAC in axial and coronal reformats in bone window 0.5mm X 0.5mm |
| Comments | |
| Images to send to PACS | Scouts, axial and coronal images in bone and soft tissue, magnified axial and coronal images of right and left side |
| | |
| | |