

BRYANRADIOLOGYASSOCIATES

2722 Osler Boulevard Bryan, TX 77802 Phone: (979) 776-8291 Fax: (979) 774-7871

Hemangioma Abdomen

Indications :

Scan Type	Helical
Scan Start\ End Location	Above diaphragm to iliac crest
Scan Details	5mm X 5mm Non contrast Arterial phase Venous phase 3 minute delay 5 minute delay 10 minute delay (if needed)
IV Contrast	70ml injected at 3cc/sec
Scan Delay	Arterial phase starts when descending aorta is highlighted or reaches 180 HU for sure start. 65 second venous delay 3 minute delay 5 minute delay 10 minute delay (if needed)
Reformat Details	Coronals on venous phase only, unless asked by radiologist to give more
Comments	IV and oral contrast always used. You don't always need to do the 10 minute delay, ask radiologist first.
Images to send to PACS	Scouts, axials from all phases, coronals from venous