

BRYAN RADIOLOGY ASSOCIATES

2722 Osler Boulevard Bryan, TX 77802 Phone: (979) 776-8291 Fax: (979) 774-7871

Carotids

Indications : stenosis, aneurysm

Scan Type	Helical
Scan Start\ End Location	1 cm below aortic arch to the circle of willis Scan from bottom to top
Scan Details	2mm X 2mm in an angio window
IV Contrast	100cc injected at 4-5cc/sec
Scan Delay	Sure start at aortic arch, trigger at 120-130 HU, or when internal carotids start blushing
Reformat Details	2mm X 2mm coronals and sagittals in angio window, volume is sent to vitrea
Comments	3D internal and external carotid images are worked up on the vitrea
Images to send to PACS	Scouts, axials, coronals, sagittals, and 3D worked up from Vitrea