

BRYANRADIOLOGYASSOCIATES

2722 Osler Boulevard Bryan, TX 77802 Phone: (979) 776-8291 Fax: (979) 774-7871

Angio Abdomen Pelvis for AAA

Indications : evaluate for aneurysm

Scan Type	Helical
Scan Start\ End Location	Above diaphragm through pelvis.
Scan Details	1mm X 1mm in angio window Noncontrast Contrast arterial phase
IV Contrast	100cc injected at 4cc/sec
Scan Delay	Use sure start, put ROI in aorta trigger at 180 HU or start the scan when the aorta is highlighted
Reformat Details	2mm X 2mm coronals and sagittals images, volume sent to Vitrea
Comments	N/A
Images to send to PACS	Scouts, axials, coronals and sagittals, 3D from vitrea