

BRYANRADIOLOGYASSOCIATES

2722 Osler Boulevard Bryan, TX 77802 Phone: (979) 776-8291 Fax: (979) 774-7871

Aortic Run-off (peripherals)

Indications :

| | |
|--------------------------|--|
| Scan Type | Helical |
| Scan Start\ End Location | Above diaphragm through feet |
| Scan Details | 2mm X 2mm in soft tissue |
| IV Contrast | 100ml injected at 4-5cc /sec |
| Scan Delay | Use sure start, put ROI in descending aorta, trigger at 180 HU or start the scan when aorta is highlighted |
| Reformat Details | 2mm X 2mm coronals and sagittal images, volume to vitrea |
| Comments | N/A |
| Images to send to PACS | Scouts, axials, coronals and sagittals 3D from Vitrea |
| | |
| | |