

BRYANRADIOLOGYASSOCIATES

2722 Osler Boulevard Bryan, TX 77802 Phone: (979) 776-8291 Fax: (979) 774-7871

Angio Chest

Indications : evaluate for aneurysm

Scan Type	Helical
Scan Start\ End Location	Above apices through top of kidneys
Scan Details	3mm X 3mm in angio window
IV Contrast	100ml injected at 4cc/sec
Scan Delay	Use sure start, put ROI in aortic arch, trigger at 180 HU or start the scan when the arch is highlighted
Reformat Details	3mm X 3mm coronals and sagittals
Comments	N/A
Images to send to PACS	Scouts, axials, coronals and sagittals