

BRYAN RADIOLOGY ASSOCIATES

2722 Osler Boulevard Bryan, TX 77802 Phone: (979) 776-8291 Fax: (979) 774-7871

Aortic Dissection

Indications :

Scan Type	Helical
Scan Start\ End Location	Above arch through bifurcation
Scan Details	3mm X 3mm in
IV Contrast	100cc injected at 4cc/sec
Scan Delay	Use sure start, put ROI in descending aorta, trigger at 180 HU or start the scan when aorta is highlighted
Reformat Details	3mm X 3mm coronals and sagittals images , volume sent to Vitrea
Comments	N/A
Images to send to PACS	Scouts, axials, coronals and sagittals, 3D images from vitrea