

BRYANRADIOLOGYASSOCIATES

2722 Osler Boulevard Bryan, TX 77802 Phone: (979) 776-8291 Fax: (979) 774-7871

Routine Pelvis

Indications : pain, constipation, mass

Scan Type	Helical
Scan Start\ End Location	Above iliac crest to greater trochanters
Scan Details	5mm X 5mm abdomen window
IV Contrast	70ml injected at 2cc/sec
Scan Delay	40 second
Reformat Details	Coronals in abdomen window
Comments	IV and oral contrast are used for this exam
Images to send to PACS	scouts, axials, coronals