BRYAN**RADIOLOGY**ASSOCIATES

2722 Osler Boulevard Bryan, TX 77802 Phone: (979) 776-8291 Fax: (979) 774-7871

Chest (thorax)

Indications :

Non Contrast nodule, Dr prefers non contrast Contrast nodule, mass, SOB, metastasis, pneumonia, abnormal CXR, cough

Scan Type	Helical
Scan Start\ End Location	Above apices through top of kidneys
Scan Details	3mm X 3mm in soft tissue and lung windows
IV Contrast	70cc injected at 2cc/sec
Scan Delay	Start scan after all contrast has been injected
Reformat Details	3mm X 3mm coronals and sagittals in soft tissue window
Comments	N/A
Images to send to PACS	Scouts, axials in soft tissue and lung window, coronals and sagittals in soft tissue window