

BRYANRADIOLOGYASSOCIATES

2722 Osler Boulevard Bryan, TX 77802 Phone: (979) 776-8291 Fax: (979) 774-7871

Chest for (low dose cancer screening)

Indications :

Scan Type	Helical
Scan Start \ End Location	Above apices through lung bases
Scan Details	2mm X 2mm in soft tissue and lung windows
IV Contrast	N/A
Scan Delay	N/A
Reformat Details	2mm X 2mm coronals
Comments	Scan parameters need to meet ACR requirements for your scanner. Lowest Dose possible is used
Images to send to PACS	Scouts, axials in soft tissue and lung window, coronals lung window