

BRYAN RADIOLOGY ASSOCIATES

2722 Osler Boulevard Bryan, TX 77802 Phone: (979) 776-8291 Fax: (979) 774-7871

Angio Chest for PE

Indications : evaluate for pulmonary embolism

Scan Type	Helical
Scan Start\ End Location	Above apices through top of kidneys
Scan Details	3mm X 3mm in soft tissue and lung windows
IV Contrast	100cc injected at 4cc/sec
Scan Delay	Use sure start, put ROI in pulmonary trunk, trigger at 180 HU or start the scan when the pulmonary trunk is highlighted.
Reformat Details	3mm X 3mm coronals and ...oblique MIP images need to be performed for the right and left pulmonary arteries
Comments	N/A
Images to send to PACS	Scouts, axials in soft tissue and lung window, coronals and oblique MIPS in soft tissue window