

BRYAN RADIOLOGY ASSOCIATES

2722 Osler Boulevard Bryan, TX 77802 Phone: (979) 776-8291 Fax: (979) 774-7871

Enterography with VoLumen (abdomen wwo pelvis wwo)

Indications : Crohn's Disease, evaluate small bowel

Scan Type	Helical
Scan Start\ End Location	Above diaphragm to greater trochanters
Scan Details	5mm X 5 mm abdomen window Non contrast 40 second delay 70 second delay
IV Contrast	70ml injected at 2cc/sec
Scan Delay	40 second delay after start of contrast 70 second delay
Reformat Details	Coronals and sagittals on 40 sec and 70 sec delay in abdomen window
Comments	Pt drinks 3 bottles of VoLumen prior to exam and needs to be scanned within 45 minutes of drinking the first bottle. IV contrast has to be used for this exam otherwise it is not diagnostic.
Images to send to PACS	Scouts, axials from all 3 phases, coronals and sagittals