

# BRYANRADIOLOGYASSOCIATES

2722 Osler Boulevard Bryan, TX 77802 Phone: (979) 776-8291 Fax: (979) 774-7871

## Bony Pelvis

**Indications :** Fracture

Scan Type	Helical
Scan Start\ End Location	Above iliac crest to greater trochanters
Scan Details	5mm X 5mm bone window
IV Contrast	N/A
Scan Delay	N/A
Reformat Details	Coronals and sagittals in bone window
Comments	N/A
Images to send to PACS	scouts, axials, coronals, sagittals